

Visit us : www.ashbme.com

E-mail : monil.bmed@gmail.com



BIO-MEDICAL INSTRUMENT'S SERVICES

AUSSIN SHUBH HOSPICA PVT. LTD.
SF-8, JBR ARCADE, SCIENCE CITY ROAD, AHMEDABAD,
GUJARAT, INDIA-380061

SALES | SERVICE | AMC-CMC | CALIBRATION | REFURBISHED

PH : +91 79 65401100, MOBILE : +91 99781 63777, +91 9925063777

MODULAR OPERATION THEATER EQUIPMENTS



Our Company Profile

We would like to introduce ourselves as a professionally managed company, pioneer in the business of door hardware and Accessories for


- HVAC
- Cold rooms
- Clean rooms
- Modular operation theater

We are a Govt. recognized leading manufacturer and exporter of PVC Corner Coving, Mortise Locks, SS Hinges, SS handles, Nylon corners, Window View ports, Rubber Gaskets, PVC gaskets, aluminum die casted hinges, etc & These Products are made of quality raw material.

We are having upbeat infrastructure and skilled staff. Apart from that, we are very cautious to maintain the quality of the products, which is controlled by top of the line professionals whose excellent performance brought laurel to the company placing it in the forefront among the top ranking manufacturers & exporters in India.

We are committed to the customer satisfaction by supplying door hardware and Accessories for clean rooms, telecom shelters, air handling units at mutually agreed terms. We strive for continuous improvement in Quality Assurance Systems, productivity and enhance the skills of all it's personnel.


We are presently exporting to various reputed buyers in Europe and Asian countries. Moreover we also have a well-managed unit, where we have manufacture our range of products.



OT CONTROL PANEL

OT Control Panel


Backed by the profound experience of the market, we are offering a commendable range of Surgeon Control Panel. The products we are offering are designed to scope with changing technology and equipment in operating environments. User friendly, high performance, smooth finishing, longer service life, ease of operating and maintaining purpose are some of its salient features. Our patrons can avail this range as per their specifications



ALT-608

LED X-RAY VIEW BOX

LED X-Ray View Box
Product code:- ALT-609



HERMATIC SEALED SLIDING DOOR




ALT-607

FIXED PENDENT

Our firm has emerged as the most remarkable name in this domain involved in offering the best quality array of Fixed Pendant that is accurately designed according to the set global parameters. Widely used in operation theatre and is available in assorted specifications according to the varied application needs of hospitals. Apart from this, offered product is strictly tested for several flaws before making their final delivery at customers' premises. Features:

- Easy to operate
- Easy installation
- Durable finishing



ALT-605

MAGNETIC WRITING LIST BORAD



ALT-610

SURGICAL SCRUB SINK

ALTOS offers scrub stations which are designed to ensure that surgeons and staff can undergo a thorough aseptic scrub, whether using the count-stroke or timed scrub methods. The scrub sink is fabricated from type 304 heavy gauge stainless steel. The stainless steel sink is polished to a satin finish.

The Scrub Sink is provided with a tiltable top for access to the water control valves, waste connections, stops and strainers. The scrub sink is fitted out with quality tap with external sensors and automatic sensor operated Betadine dispensers. The unit is designed to be fixed on wall using metal brackets or floor mounted.

In TIMER model User can set the scrub time from 1 to 10 minutes 6 hours continuous battery backup. A thermostatically controlled mixing valve with anti-scald system automatically maintains water temperature.



ALT-604



PRESSURE RELIEF DAMPER

We at, Altos engineers are engaged in manufacturing and exporting Pressure Relief Damper. It is manufactured in compliance with international standards. This product is featured with unique design, which is capable of preventing the operation theaters from contamination of air dirty area. These are provided with suitably sized air pressure relief damper placed, to enable differential room pressure to be maintained the room infection-free and sterile.



ALT-613

LAMINAR AIRFLOW

Laminar Airflow is designed to move particle free air over the aseptic operation field in one direction. It can be designed to flow vertically or horizontally and is usually combined with high efficiency particulate air (HEPA) filters. HEPA filters remove particles >0.3 micron in diameter with an efficiency of 99.97%.



ALT-606

Visit us : www.ashbme.com E-mail : monil_bmed@gmail.com

ASH REPAIRERS

Aussin

BIO-MEDICAL INSTRUMENT'S SERVICES

AUSSIN SHUBH HOSPICA PVT. LTD.
SF-8, JBR ARCADE, SCIENCE CITY ROAD, AHMEDABAD,
GUJARAT, INDIA-380061

SALES | SERVICE | AMC-CMC | CALIBRATION | REBURNISHED

PH : +91 79 65401100, MOBILE : +91 99781 63777, +91 9925063777

HOSPITAL SURGICAL PENDANT



ALT-603

Anesthesia Pendant

We hold great expertise in providing the best array of Anesthesia Pendant. Offered product is used in hospitals for giving anesthesia to the patients. Perfectly designed by making use of the best quality raw material procured from the certified vendors of the market, this assures long life. In addition to this, to meet diverse needs of our customers, we offer this product in assorted specifications. Features:

Robustness

Rust resistant

Longer service life



PASS BOX HATCH BOX



ALT-614

Pass Box Hatch box

We are amongst the reputed enterprises, engaged in providing the best quality array of Pass Box hatch box. Offered product is precisely developed as per the guidelines defined by our talented experts; offered product is accurately manufactured making use of optimum grade raw material that is procured from some of the certified vendors of the market.

Size:- coustmise

Material – SS-304

LED PERIPHERAL LIGHT



ALT-611

LED Peripheral Light

We offer 2nd Generation LED Surgical light, innovative thought to build up operating room with powerfull illumination light for better surgical. Our LED Surgical light provides cold, shadowless and brighter light area for long surgeries, without any interruption to surgeons. LEDs have a technical advantage over the traditional halogen lamp, which produce heat, limited life cycle of bulb. Thanks to LEDs technique, now surgeon can work in low heat and great life cycle of LEDs working hours, with low power consumption & natural white day light.



ALTOS ENGINEERS PVT.LTD.

Works: Plot No 758, Street No. 13, Kadipur Ind. Area, Gurgaon Haryana (India)

Mob.: +91-8750050499, 9711360601

E-mail : altosindia@gmail.com, www.altosengineers.com